

# BOLSACREPE



Excellent quality crepe paper in green colour. Soft and malleable material, to easy fold it without any marks



For single use and for double wrapping

## GREEN CREPED PAPER

Green creped paper sheets of 60 g for wrapping and sterilization of medical products. It is specially used in hospitals, medical centres, veterinarian centres, CSSD, etc.

It is also used as tray protector and for double wrapping. Manufactured from 100% cellulose and conforms to ISO 11607-1 and EN 868-2. Basically for wrapping medical devices to be sterilized.

It is considered a single use product and can be sterilized by Steam, Ethylene Oxide, Formaldehyde or Irradiation.

## CHARACTERISTICS

- Soft and malleable
- Excellent bacterial barrier
- Liquid repellent
- Wide range of sizes

**bolsaplast**medical

This product is available in the following sizes:

Code	Sizes		U/Box	Box/cm	Box/Kg
MCREPV50	50 x 50 cm	19.69" x 19.69"	500	51 x 26 x7	7,6
MCREPV60	60 x 60 cm	23.62" x 23.62"	500	61 x 38 x14	10,6
MCREPV75	75 x 75 cm	29.53" x 29.53"	250	78 x 29 x 12	8,3
MCREPV90	90 x 90 cm	35.43" x 35.43"	250	91 x 50 x 10	12
MCREPV100	100 x 100 cm	39.37" x 39.37"	250	101 x 54 x 10	15
MCREPV120	120 x 120 cm	47.24" x 47.24"	100	121 x 33 x 10	8,6